



CENTRIC MEDICAL™

*Revolutionizing
FOOT AND ANKLE*

Product Description

The **ROGUE™ Hammertoe Correction System** offers a dual threaded design and is inserted between the proximal and middle phalanges, so that the opposing threads fixate on the phalangeal canal of the toe and compress the joint.

Key Features & Benefits **ROGUE & ROGUE+**

- » Compared to current standard of care, may lead to lower infection rate and less pain due to a fully implanted device
- » Dual threaded distal and proximal components for bone purchase and compression across the joint
- » Follow up surgery is not required

ROGUE+ Additional Features & Benefits

- » Less implant extruding into the PIP Joint provides easier reduction, allowing for less bone removal while maintaining original toe length
- » Reverse thread design simplifies the technique to consistently place the implant across the PIP Joint, optimizing bone purchase in the proximal and middle phalanx
- » Streamlined instrumentation for reduced OR time



ROGUE

Hammertoe Correction System



ROGUE+

Hammertoe Correction System



**CENTRIC
MEDICAL™**

ESTABLISHING THE FOUNDATION


Centric Medical is dedicated to improving the quality of life for patients with distal extremity symptomatology, increasing procedural efficiency and efficacy through design, uncompromising quality standards, and the most technologically advanced manufacturing platforms.



ROGUE
Hammertoe Correction System

ROGUE™ + Hammertoe Correction System

- 114-188 Centric Medical Tray
- 154-451 ROGUE Implant
- 154-455 ROGUE Step Drill
- 154-457 Small ROGUE Driver, 2.0mm
- 154-475 ROGUE Handle Body
- 154-476 ROGUE Handle Turn Knob



ROGUE™ Hammertoe Correction System

- 114-188 Centric Medical Tray
- 154-450 ROGUE Implant
- 154-455 ROGUE Step Drill
- 154-456 Large ROGUE Driver, 2.9mm
- 154-457 Small ROGUE Driver, 2.0mm
- 154-475 ROGUE Handle Body
- 154-476 ROGUE Handle Turn Knob
- 154-477 ROGUE Counter Push

